

# Specifications Noorlander Largo

## Keyboards

5 octaves: C-c4, 61 keys, with tracker simulation	•
3, plastic keys with steel key frame	•
3, wooden key centres with plastic covering	Option
3, wooden key centres with ebony or cherry naturals and oak, cherry or ebony sharps	Option

## Pedal

30-tone: C-f1, straight, oak keys and sharps	•
30-tone: C-f1, concave, oak keys and sharps	Option
Wenge (dark) sharps	Option

## Swell pedals

2 pieces	•
Solid oak	•

## Touchscreens

Number of touchscreens	2
Size touchscreens in inch	10

## Furniture finish

Oak wood	•
Oil (your colour)	•
Transparent varnish	Option
Coloured stain	Option
Painted	Option
Other wood type	Option

## LED lighting

Pedal lighting	•
Music rack lighting	•
Dimmed lighting	•

## Thumb pistons

14 pieces (S, 0, 1, 2, 3, 4, <, >, 5, 6, 7, 8, 9, C)	Option
--	--------

## Foot pistons

2 pieces, gold/silver colour	Option
------------------------------	--------

## Bench model

'Standard' bench	•
Storage space	•
Adjustable height	•
'German' bench	Option
'French' bench	Option
Leather covered seating (choice of colours)	Option

## Menu functions

Combinations (per sample set)	260
Selection temperaments	15
Set	•
Cancel	•
Selection sample set	•

## Audio

Number of audio channels	6
Number of woofers (middle tones)	2
Number of tweeter (high tones)	2
Subwoofer (low tones)	2

## External connections

USB 2.0	•
Headphone connection	•

## Console dimensions

Height	123 cm
Width	133 cm
Depth including pedal (ca. 25 cm pushed in)	110 cm
Depth without pedal and keyboards	65,5 cm

Specifications are subject to change



Largo



## Why a Largo?

An organ should arouse your senses. The sound is very important but so is the organ console itself. Noorlander Organs leads in this field. Our organs have a very realistic (pipe) organ sound, thanks to extremely advanced Sweeling technology developed by Noorlander Organs.

For optimal use of Sweeling we have developed special hardware and software.

Our organs have a fantastic sound reproduction system aided by a Noorlander Organ Manager. This software automatically regulates the power, audio and midi. The finishing touches can be done to match your other furniture. All operating controls are carefully disguised in the organ console.

We can offer you personalized advice and a premium guarantee. Do you need help? Through internet we can offer you the help you need.



### About Sweeling

#### Technology

Sweeling is the technology we have developed to play virtual pipe organs. The biggest advantages of Sweeling are the infinitely adjustable reverberation, the super-fast loading time, the simplicity of operation, the very realistic wind model and the great sound realism. The organs are recorded pipe by pipe. Details such as the keyboard and stop action are also included. The sound collection of the recordings of these pipes is called a **sample set**.

#### Organs

Famous organs have since been recorded for Sweeling, such as the organs of the Grote Kerk in Dordrecht (NL). The sound fidelity can best be compared to a true-to-life CD recording. Sweeling and its technique appeal to the imagination, the possibilities are almost limitless! Furthermore, numerous sample sets that are already on the market are suitable for playing in Sweeling.



## Largo

### Sound

- High-end digital audio system
- A much higher quality sound than other electronic organs
- The sound of a real pipe organ that is bound to increase your playing pleasure

### Usage

- Very easy to use
- Easy and quick stop changes with the use of the touchscreens
- You can program the stop combination pistons yourself
- You don't need to be a computer expert. The menu is very simple to use.



Largo:  
"A compact classic..."



### Design

- Classic design with lots of solid oak wood
- Compact console
- Very high quality finishing and beautiful panelling
- Standard LED lighting in music rack and above pedalboard

### Possibilities / additions

- Easy to expand with other (pipe organ) sample sets
- Other sample sets and temperaments easily uploaded
- Easy to add a surround sound system